



Mapping for Chronic 64 Channel Probe

ASSY-156 ECoG-64A



Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



48	50	1	31
47	17	15	49
18	29	3	16
32	45	51	2
46	20	14	52
19	30	4	13
42	43	53	56
44	22	12	54
28	39	57	5
40	41	55	58
23	28	6	9
36	37	59	62
38	24	10	60
25	33	63	7
34	35	61	64
21	26	8	11