



Mapping for Chronic 64 Channel Probe

ASSY-156 H1



Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



Shank A				Shank B			
26	33	37	41	12	4	51	15
35	24	28	39	53	14	3	1
27	40	23	36	31	49	16	2
38	25	34	21	52	13	56	54
22	43	30	20	8	61	63	10
45	29	17	50	59	6	55	57
48	47	18	32	5	58	9	62
46	19	42	44	60	7	64	11