



Mapping for Chronic 64

Channel Probe

ASSY-156 M1 (June 2022 Onwards)



Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



26		36	
35	41	38	43
33	39	25	30
24	27	34	20
37	40	21	45
28	23	22	29
48	17	47	50
46	18	19	32
12	42	53	44
51	4	3	14
31	15	49	1
52	16	13	2
8	56	61	54
59	63	6	10
5	55	58	57
60	9	7	62
	64		11

