



Mapping for Chronic 64

Channel Probe

ASSY-156 M2 (June 2022 onwards)



Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



26	12
35	53
33	4
24	14
37	51
28	3
41	15
39	1
27	31
40	49
23	16
36	2
38	52
25	13
34	56
21	54
22	8
43	61
30	63
20	10
45	59
29	6
17	55
50	57
48	5
47	58
18	9
32	62
46	60
19	7
42	64
44	11

Upper-half

Lower-half

