




Mapping for Y Chronic 64 Channel Probe

ASSY-158 E-1 & E-2

Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



Shank A

Shanks B - D

61	48	18	3
60	45	19	6
59	34	32	5
58	47	17	8
62	46	20	4
57	43	21	7
56	44	22	10
55	33	31	9
63	41	23	1
54	42	24	12
53	39	25	11
52	36	30	14
49	40	26	15
64	37	27	2
51	38	28	13
50	35	29	16