




Mapping for Chronic 64-Channel Probe

Molex Connector ASSY-236 E-1 & E-2

Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).

TOP

BOTTOM



Shank A	Shank B	Shank C	Shank D
64	46	34	49
63	45	39	50
48	30	18	19
15	11	6	2
16	58	54	1
62	44	36	51
61	27	23	35
31	28	20	37
13	10	7	4
14	57	55	3
60	42	38	52
47	43	41	33
29	26	22	21
32	9	8	17
12	24	40	5
59	25	56	53

Shank A

Electrodes
this side