



# Mapping for Chronic 64-Channel Probe

## Molex Connector ASSY-236 H10

Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).

**TOP**

**BOTTOM**



	46	
16		44
	63	
48		11
	30	
64		45
	58	
61		28
	27	
13		57
	10	
60		43
	42	
29		9
	15	
12		25
	62	
59		24
	31	
32		26
	14	
	47	

**Shank A**

	49	
54		51
	39	
18		2
	19	
34		50
	1	
23		37
	35	
7		3
	4	
38		33
	52	
22		17
	6	
40		53
	36	
56		5
	20	
8		21
	55	
	41	

**Shank B**

Elect this