

Mapping for Digital 128-Channel Probe

SPI connector

ASSY-350 H14-1 / H14-2

(was ASSY-350 H11-1 / H11-2 April - June 2024)

Electrodes are located on the side of the probe as indicated by the arrow below - use Shank A orientation as a reference to construct the channel-map with respect to your brain target. Electrode numbers are 0-indexed so as to correspond to the Intan / Open Ephys data-acquisition software channels.

Red Wire and White Wire correspond to indicated REF and GND inputs to the headstage amplifier.

