



Mapping for Acute 32-Channel Probe

ASSY-37 H10b



Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



Electrodes
this side

	17	
8		19
	15	
12		23
	21	
16		18
	22	
6		20
	25	
2		24
	27	
3		26
	29	
5		31
	10	
9		32
	14	
11		30
	13	
7		28
	4	
	1	

