




Mapping for Acute 64-Channel Probe

ASSY-77 E-1 & E-2

White Wire
GND

Red Wire
REF

Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



	Shank A	Shank B	Shank C	Shank D
21	1	64	44	
22	2	63	43	
23	3	62	42	
28	4	61	37	
24	5	60	41	
32	6	59	33	
30	7	58	35	
26	8	57	39	
29	9	56	36	
25	10	55	40	
16	11	54	49	
17	12	53	48	
18	13	52	47	
19	14	51	46	
20	15	50	45	
27	31	34	38	