



Mapping for Acute 64-Channel Probe

ASSY-77 H1

White Wire
GND

Red Wire
REF

Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



Shank A

Shank B

21	24	29	18
23	30	16	20
27	19	17	25
26	32	28	22
1	3	5	7
9	11	13	15
31	14	12	10
8	6	4	2

64	60	56	52
62	58	54	50
34	51	53	55
57	59	61	63
44	42	41	35
36	49	47	45
38	46	48	40
39	33	37	43

