


# Mapping for Acute 64-Channel Probe

ASSY-77 H6

White Wire  
GND

Red Wire  
REF

Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



## Shank A

## Shank B

	17
21	25
23	26
24	32
30	28
29	22
16	1
18	3
20	5
27	7
19	9
11	13
15	31
14	12
10	8
6	4
2	

	49
53	57
55	58
56	64
62	60
61	54
48	33
50	35
52	37
59	39
51	41
43	45
47	63
46	44
42	40
38	36
34	

