


Mapping for Acute 64-Channel Probe

ASSY-77 M1 (Discontinued May 2022)



White Wire
GND

Red Wire
REF

Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



21	25	
	18	3
23	26	
	20	5
24	32	
	27	7
30	28	
	19	9
29	22	
	17	11
16	1	15
	13	14
31	12	10
	8	6
8	4	2
	63	61
	59	57
	55	53
	51	34
	50	52
	54	56
	58	60
	62	64
	43	37
	33	39
	40	48
	46	38
	45	47
	49	36
	35	41
	42	44



Electrodes
this side