


# Mapping for Acute 64-Channel Probe

ASSY-77 M1 (June 2022 Onwards)



**White Wire**  
GND

**Red Wire**  
REF

Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



21	25	
	18	3
23	26	
	20	5
24	32	
	27	7
30	28	
	19	9
29	22	
	17	11
16	1	
	13	15
31	14	
	12	10
8	6	
	4	2
64	62	
	60	58
56	54	
	52	50
34	51	
	53	55
57	59	
	61	63
44	42	
	41	35
36	49	
	47	45
38	46	
	48	40
39	33	
	37	43



Electrodes  
this side