


Mapping for Acute 64-Channel Probe

ASSY-77 M2 (Discontinued May 2022)

White Wire
GND

Red Wire
REF

Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



- 21
 - 23
 - 24
 - 30
 - 29
 - 16
 - 18
 - 20
 - 27
 - 19
 - 17
 - 25
 - 26
 - 32
 - 28
 - 22
 - 1
 - 3
 - 5
 - 7
 - 9
 - 11
 - 13
 - 15
 - 31
 - 14
 - 12
 - 10
 - 8
 - 6
 - 4
 - 2
- Upper-half**
- 63
 - 61
 - 59
 - 57
 - 55
 - 53
 - 51
 - 34
 - 50
 - 52
 - 54
 - 56
 - 58
 - 60
 - 62
 - 64
 - 43
 - 37
 - 33
 - 39
 - 40
 - 48
 - 46
 - 38
 - 45
 - 47
 - 49
 - 36
 - 35
 - 41
 - 42
 - 44
- Lower-half**

