


Mapping for Acute 64-Channel Probe

ASSY-77 M2 (June 2022 Onwards)

White Wire
GND

Red Wire
REF

Red Wire and White Wire correspond to indicated REF and GND pins on the probe connector - verify how these connect to your headstage pre-amplifier (1 wire may be sufficient if you prefer GND / REF shorted-together in your headstage).



21	64
23	62
24	60
30	58
29	56
16	54
18	52
20	50
27	34
19	51
17	53
25	55
26	57
32	59
28	61
22	63
1	44
3	42
5	41
7	35
9	36
11	49
13	47
15	45
31	38
14	46
12	48
10	40
8	39
6	33
4	37
2	43

Upper-half

Lower-half

