


Mapping for Deep Array 64-Channel Probe

SPI connector

DA 64-2

Electrode numbers are 0-indexed so as to correspond to the Intan / Open Ephys data-acquisition software channels.



Upper-half

	32
31	
	33
30	
	34
29	
	35
28	
	36
27	
	37
26	
	38
25	
	39
24	
	40
23	
	41
22	
	42
21	
	43
20	
	44
19	
	45
18	
	46
17	
	47
16	



Lower-half

	48
15	
	49
14	
	50
13	
	51
12	
	52
11	
	53
10	
	54
9	
	55
8	
	56
7	
	57
6	
	58
5	
	59
4	
	60
3	
	61
2	
	62
1	
	63
0	

