


Mapping for Mini-Amp-64 V1 Headstage

Molex Connector

Compatible with ASSY-236



Connector J1
Probe ("BOTTOM")

INTAN	Probe	INTAN	Probe
-		-	
63	49	55	33
62	50	53	34
61	51	54	35
60	52	52	36
59	53	51	37
58	54	50	38
57	55	49	39
56	56	48	40
7	57	47	41
6	58	15	42
5	59	13	43
4	60	12	44
3	61	11	45
2	62	9	46
1	63	10	47
0	64	8	48

Connector J2
Probe ("TOP")

INTAN	Probe	INTAN	Probe
REF	REF	GND	GND
46	17	41	1
45	18	39	2
44	19	38	3
43	20	37	4
42	21	35	5
40	22	34	6
36	23	33	7
31	24	32	8
27	25	29	9
23	26	30	10
21	27	28	11
18	28	26	12
19	29	25	13
17	30	24	14
16	31	22	15
14	32	20	16